

Easy Fluor-Western Blocking Buffer

Overview of Easy Fluor-Western Blocking Buffer:

Nepenthe Easy Fluor-Western Blocking Buffer is specifically optimized to enhance specific antibody-antigen interactions for fluorescent Western Blots. This all-in-one blocking solution is designed to improve sensitivity and decrease overall background. Non-specific binding caused by low quality antibodies is reduced while signal from the specific antibody-antigen complex is stabilized and enhanced.

Easy Fluor-Western Blocking Buffer is provided as a convenient ready-to-use 1x solution intended to directly replace other commonly used blocking buffers for fluorescent Western Blotting.

- **ENHANCED SIGNAL** – Enhances Western blotting signals
- **DECREASED BACKGROUND** – Decreased non-specific bands and lowered background
- **SENSITIVITY** – Get better signals for low protein sample concentrations
- **CONVENIENCE** – Save preparation and optimization time with pre-made buffer

For Orders:

Catalog Number	Product	Size
NB0503729A2	Easy Fluor-Western Blocking Buffer	1 L

